

Reflex Sheet– P1/2






Configure a new RESCUE event

1. Create the event

- Events
- Organizations
- Users

Home > Events > List

Events » List events

stl	All	All	All		
Type	Wording	Start	End	Organization	Actions
Trail running	Formation STL	5/14/19, 10:56 AM	11/16/22, 7:24 PM	logicoss Démonstration	   
 Create a new Event from the list of events					Item 1 to 1 of 1

- Events
- Organizations
- Users

Home > Events > Create

Events » Create an event

Organization *

Wording *

Type

Version *

The event can be editable with API ☐

Display the severity of pathologies ☐

Event Id API After

Number of expected spectators

Only inform the spectators here. Attendees will also be counted since the import of attendees listings.

Logo PNG or JPEG picture, from 25x40px to 300x100px, max 300 Ko.

Time zone *

Start *

End *

Latitude *

Longitude *

Map style *



Save

Back to the list Save

2. Creation of a second aid station or more

Notice: the first aid station is created by default when the event is created!

- Events
- Organizations
- Users

Home > Events > Formation logicoss Rescue > First aid stations > Create

First aid stations » Create first aid station(s)

Formation logicoss Rescue

Event Settings Medical Form menu.Interfaces Evacuation centers Files Import participants KMLs Intervention Assets First aid stations Landmarks SMS Content Users Map edition View

Wording *

Latitude *

Longitude *

Family

Phone


Note(s)

beds

Wording *

Latitude *

Longitude *



3. Creation/Assignment of users

Case 1: The user already exists AND has already been associated with the rescue service’s organization

Home > Events > Mask > Users > List

Users > List users

Mask

Event Settings Medical Form menu.interfaces Evacuation centers Files Import participants KMLs Intervention Assets First aid stations Landmarks SMS Content Users Map edition View

logicoss users

Common name	Common first name	Mobile	Type	Family	First aid station	Profile	Last connection
No element							

New user

E-mail address	Common name	Common first name	Mobile	Type	Family	First aid station	Profile
No element							

Edit

1 Select permanent users from the list

2 Restrict the access to a single aid station in the list

3 Assign a profile "01. RESCUE Event Staff"

4 Add other users or save

Case 2: The user doesn't yet have an account or is not associated yet on the organization

logicoss users

Username * Type * Family First aid station * Profile *

George SWED Administrative staff ... Poste de secours Principal No right

New user

E-mail address *	Common name *	Common first name *	Mobile	Language	Type	Family	First aid station *
g.swed@logicoss.com	SWED	George	(+33)66134325	Anglais	Administrat...	...	Poste de secours 1

01. RESCUE Event Staff

02. RESCUE Event Admin

03. RESCUE Event Omniscient

08. RESCUE Families Information Point (ex Info Families)

No right

Back to the list Save

5 Each new user receives an e-mail allowing him to set his password

1 Fill in the surname, first name, mobile number in the "New User" form & language

2 Possibly restrict access to a single aid station

3 Grant a "01. RESCUE Event Staff"

4 Save

OR

CREATE YOUR PASSWORD

A registration in your name has been made on the platform Logicoss. Therefore please click on the following link to create your password. You can then connect with the following login : utmbtraining

https://logicoss.net/reset-password/reset/uYOHJcy892aPa6wLFcEuCqkyBNXTnDQ6cITfTC/1


Logicoss Digital Healthcare, 109 Boulevard de l'Europe 69310 PIERRE BENITE, FRANCE








+33 (0)4.72.30.01.03











contact@logicoss.com

User Manual - P1/1

Download the PDF Report of the event

★ Demonstration WA - Marathon de Rouen
9/25/19, 12:51 AM - 10/14/19, 12:51 AM

13/51/417/1815/370/1Gaspard LEBEL

InterventionsFirst aid stationsAssetsMarshalsBedsAttendeesMapOccurrence booksReportsSettings

Reports » Report management

Create a report Generate

Report profiles :

Parts available for report generation :



- ☒ Synthetic report
 - ☐ Description of the first response organization deployed
 - ☐ Specific comments
 - ☒ Overall activity report
 - ☒ Activity report sorted by reason
- ☒ Real-time report
 - ☐ Primary and secondary flows
 - ☒ Solicitation of rescue assets outside the organization
 - ☒ Solicitation of the evacuation centers
 - ☐ Main reasons of consultations
- ☒ Interventions
 - ☒ Interventions status
 - ☒ Reasons of Interventions' completion
- ☐ First aid station consultations
 - ☐ Distribution of consultations in first aid stations
 - ☐ Reasons for consultation in first aid stations
 - ☐ Completion reasons in first aid stations
- ☒ Accidentology
 - ☐ Accidents Map
 - ☒ Accidents distribution per landmarks
 - ☒ World Health Organization (WHO) classification
 - ☐ PPR, RTHR, TTHR, ITR
- ☒ Pathologies
 - ☒ Number of pathologies filled
 - ☒ Location of the pathologies filled
- ☒ Evacuations
 - ☐ Patients identity
 - ☒ Service

Period to study :
09/25/2019 12:51 AM - 10/14/2019 12:51 AM

Periodicity of division of tracking graphs :
None: only one graph for the entire period

Specific comments :

A T B I S U

List reports								
Timestamp	Creator	Start	End	Periodicity	Requested parts	Status	Actions	
10/2/19, 5:15 PM	Gaspard LEBEL	9/25/19, 12:51 AM	10/7/19, 12:51 AM	None: only one graph for the entire period	Synthetic report Real-time report Accidentology Pathologies Evacuations	Published		

Download the report as soon
as it has been generated!

User Manual - P1/1

Download an Health Folder

First aid stations

1 **+** **Create**

Consultations list

ON SITE	1:03 PM	55	
19	11/24/21		
ON SITE	7:12 PM	652	Carola DAHLMANN DRUG
18	3/3/20		
ORIENTATION	12:33 PM	698	Doro HEGGENBERGER
5	9/28/19		

Items 1 to 3 of 3 **View all**

2 **Open the health folder**

19 ON SITE

Medical advise or Orientation **3** **Declare consultation completion** **Save**

Consultation of 11/24/21 at 1:03 PM

Mariella RICKENBACHER (29 ans)

3 **Click on the identity**

Circumstances **Check-up** **Medical history** **Files** **Civil status** **Orientation** **Completion** **Care path** **1**

Notes

Creator *	Timestamp *	Origin *	Message *

Attendees > View attendee "Mariella RICKENBACHER"

Civil status

Identity

1

Birth name

Common name Rickenbacher

Common first name Mariella

Date of birth 12/4/91

Age 29 year(s)

Nationality BEL (Belgium)

Geographic coordinates

Usual address Way 3632 Villa 2221 Province :Oman
120 OM

Contact details

Context

Type Attendee

Bib 55

Various

Profession

Insurance

Note(s)

5 **Open the PDF using Adobe for example**

6 **Save**

7 **Consult the Final Health Folder**

4 **Launch PDF printing**

Print

logicoss

03/06/2017

General info

16:45 Intervention (Noct - Blessure / Traumatisme)

17:03 Transfert vers PMA Antioche par véhicule dédié

17:23 Soins PMA Antioche (Noct - Blessure / Traumatisme)

17:28 Demande de régulation

17:29 Recusation vers Lariboisière par AMU 2

17:31 Recusé vers Lariboisière par AMU 2

Identité

Genre Femme

Nom de naissance

Nom courant RICKENBACHER

Prénoms courant Mariella

Date de naissance 12/04/1991

Age 29

Lieu de naissance

Contexte

Type Particulier

Dossier 18

Divers

Profession

Assurance Non

Notes

Coordonnées géographiques

Adresse complète temporaire

Adresse complète habituelle

Coordonnées numériques

Adresse email

Téléphone fixe

Mobilité

Contacts d'urgence

Contact d'urgence

Téléphone

Médecin traitant

Téléphone du médecin

logicoss

Blans

Créateur	Modifié	Origine	Message
DRM/218	03/06/2017 16:45	Prod	Blanc
Requet 1898	03/06/2017 17:03	Logiciel	Blanc sauvegardé. Paramètre corrigé (218).

Prescriptions

Créateur	Modifié	Origine	Molécule(s) et posologie(s)
Requet 1898	03/06/2017 17:03	Méa 1	Paracétamol 500mg

Constantes

Créateur	Modifié	Circulatoire	Respiratoire	Neurologique	Métabolique
Requet 1898	03/06/2017 17:03	17 128	20	10	1.25

Pathologies

Antidépresseurs

Créateur	Modifié	Origine	Message
Requet 1898	03/06/2017 17:03	Prod	Blanc

Traitement chronique

Créateur	Modifié	Origine	Message
Requet 1898	03/06/2017 17:03	Prod	Blanc